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
The screenshot shows the DentalCAM V7.05.00 software interface. The top bar includes the 'denta^{swiss}®' logo and the text 'dentally different'. The main interface is divided into several sections:

- Top Bar:** Contains navigation icons (home, search, etc.) and the software title 'DentalCAM V7.05.00'.
- Left Panel:** Contains three main sections:
 - 排放流程 (Emission Process):** A row of icons representing different dental procedures.
 - 新增料塊 (Add Blank):** A form for creating a new blank with fields for 'Blank type' (Zirconia -D), 'Dimensions' (Ø 98.5), 'Height' (10.0), 'Scaling' (1.250), 'Lot number', 'Comment', 'No. of colored layers' (1), and 'Tool type' (Milling tools).
 - 紀錄篩選 (Record Filter):** A filter section with fields for 'Blank type', 'Height', 'ID', 'Lot number', 'Comment', and a 'Calculable jobs' checkbox.
- Blank archive:** A central area displaying a list of blanks. The first entry is '1 Wax -D' with dimensions 'Ø 98.5' and '10.0', and a '1.000' value below it. Action icons (refresh, search, delete) are visible on the right.
- Right Panel:** A dark grey area with the label '料塊紀錄 (Blank Record)'.

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The screenshot shows the DentalCAM V7.05.00 interface. The 'New blank' panel is active, showing fields for Blank type (PMMA -D), Dimensions (Ø 98.5), Height (16.0), and Scaling (1.000). A 'Blank archive' table lists existing blanks. Callouts 1-6 describe the workflow: 1. Material selection (Zirconia, translucent -D), 2. Dimension selection (16.0), 3. Scaling (1.000), 4. Naming the blank, 5. Confirming the new blank, and 6. Calling in the blank. A 'Scaling' callout box shows 'Scaling 1.000' and '收縮比例'.

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The screenshot shows the DentalCAM V7.05.00 interface. The steps are as follows:

- 1. 檔案位置**: The file location is set to 'g:\Coping' in the file browser.
- 2. 選擇檔案**: A list of files is shown, with 'full anatomical 1.stl' selected.
- 3. 剖面檢視**: A 3D model of a crown is shown in the preview window.
- 4. 類型選擇**: The 'Fully anatomical crown/bridge' option is selected in the 'Crown/bridge (less detailed)' menu.
- 5. 精緻等級**: The 'fully anatomical crown/bridge' option is selected, indicating the highest level of detail.
- 6. 拋入待車**: The 'Stereolithografie (*.stl)' filter is selected.
- 7. 待車總覽**: The '3Shape (*.cam)' filter is selected.
- 8. 確定車削**: The 'Milling' icon is clicked to start the process.

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