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
The screenshot shows the DentalCAM V7.05.00 software interface. The top bar includes the 'denta^{swiss}®' logo and the text 'dentally different'. The main interface is divided into several sections:

- 排放流程 (Emission Process):** Located at the top left, it contains a row of icons for different dental procedures.
- 新增料塊 (Add Blank):** A panel on the left side for creating new blanks. It includes fields for 'Blank type' (Zirconia -D), 'Dimensions' (Ø 98.5), 'Height' (10.0), 'Scaling' (1.250), 'Lot number', 'Comment', 'No. of colored layers' (1), and 'Tool type' (Milling tools).
- 紀錄篩選 (Record Filtering):** A panel below the 'Add Blank' section with filters for 'Blank type', 'Height', 'ID', 'Lot number', 'Comment', and 'Calculable jobs'.
- 料塊紀錄 (Blank Record):** A large central panel titled 'Blank archive' showing a list of records. The first record is: '1 Wax -D Ø 98.5 10.0 1.000'. It includes icons for search, refresh, and delete.

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The screenshot shows the DentalCAM software interface with the following components and callouts:

- 1. 材料選擇 (Material Selection):** A dropdown menu on the left showing 'Zirconia -D', 'Zirconia, translucent -D' (highlighted), 'PMMA -D', and 'Wax -D'.
- 2. 尺寸選擇 (Dimension Selection):** A vertical list of dimensions from 10.0 to 30.0, with 16.0 highlighted.
- 3. Scaling:** A callout box showing 'Scaling' set to '1.000' with the label '收縮比例' (Scaling Ratio).
- 4. 料塊取名 (Blank Naming):** A callout pointing to the 'Comment' field in the 'New blank' section.
- 5. 確認新增 (Confirm New):** A callout pointing to the '+' icon in the 'Tool type' section.
- 6. 叫入料塊 (Call in Blank):** A callout pointing to the 'Add' icon (person with arrow) in the 'Blank archive' section.
- 刪除料塊 (Delete Blank):** A callout pointing to the trash icon in the 'Blank archive' section.

The 'New blank' section contains the following fields:

Blank type	PMMA -D
Dimensions	Ø 98.5
Height	16.0
Scaling	1.000
Lot number	
Comment	
No. of colored layers	1
Tool type	Milling tools

The 'Blank archive' section contains the following table:

ID	Blank type	Dimensions	Height	Scaling	Actions
3	PMMA -D	Ø 98.5	16.0	1.000	[Add] [Delete]
2	Zirconia, translucent -D	Ø 98.5	10.0	1.250	[Add] [Delete]
1	PMMA -D	Ø 98.5	10.0	1.000	[Add] [Delete]
1	Wax -D	Ø 98.5	10.0	1.000	[Add] [Delete]

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The screenshot shows the DentalCAM V7.05.00 interface. The top toolbar has a blue box around the 'Import' icon (1). The file explorer on the left shows a folder named 'Copings' (1). The file list on the right shows several .stl files, with 'full anatomical 1.stl' selected (2). The 'Filter' dropdown is set to 'Stereolithografie (*.stl)' (6). The main 3D view shows a blue dental model (3). A dropdown menu is open, showing various dental part types, with 'Fully anatomical crown/bridge' selected (4). The 'Finishing level' dropdown is set to '3 stars' (5). The '3Shape (*.cam)' button is highlighted in the bottom right (8).

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